


Q factor

By: Kyle Woo

A dark blue diagonal gradient bar that starts from the bottom left and extends towards the top right, covering the lower half of the slide.

What is Q factor?



$$Q = \left(\frac{f_0}{\Delta f} \right)$$

f_0 = resonant frequency

Δf = bandwidth

Definition: Measures energy storage vs. energy loss in resonant systems.

Meaning of High Q:

- Sharp, narrow peak in frequency response.
- Indicates low energy loss (good for MRI coils)

-Meaning of Low Q:

- Broad peak in frequency response.
- Indicates higher energy dissipation

Significance of Q: Determines system efficiency and frequency selectivity.


My Q factor testing

Objective:


- To measure the Q-factor of future coils.
- To compare the results from the NanoVNA with those obtained from a bench-top VNA for accuracy and consistency.
- To establish a reliable method for testing coil decoupling and quality.

How I did it:

Placed a decoupled loop above a coil with an approximated Q factor of 90-100



Another example for decoupling





- The next experiment aims to find the minimum overlapping distance between loops to minimize coupling (S_{21}).



- This ensures the loops are effectively decoupled, allowing independent operation and optimal performance of the coils.

- We are trying to highlight the importance of decoupling in applications like MRI coil design, where minimal interference is critical for image quality.



5cm





4cm



3.5cm



0cm



Results

- Tested loop spacing at 5 cm, 4 cm, 3.5 cm, and 0 cm to evaluate coupling effects on Q factor.
- 3.5 cm spacing produced the highest Q factor efficiency, indicating optimal coil decoupling.
- At 5 cm and 4 cm, coupling was still present, reducing Q factor performance.
- At 0 cm overlap, coils showed strong coupling and significant signal interference.
- Results confirm that 3.5 cm is the optimal overlap distance for minimizing interference and maximizing MRI coil efficiency.
- Findings will guide future coil decoupling designs to improve image clarity and reduce energy loss.